

PRODUCT BROCHURE

OPTICS OBB-A1382

Date: 19.04.2025

ENGLISH



Trinocular lens barrel

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

Microscope tube OBB-A1382

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1382
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761176119
Customs tariff number	90119000

Technical data - Primary

Tube	Trinocular
------	------------

Technical data - Construction

Tube type	Butterfly
Tube angle	30°
Tube rotation	yes
Eye distance (max)	75 mm
Eye distance (min)	50 mm
Adjustable dioptre	both-sided

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×185×165 mm
------------------------------	----------------

Gross weight	2 kg
Net weight	1.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	600 g
Delivery time	1 d

FUNCTIONS

Standard

