

PRODUCT BROCHURE

KERN 317-01-2020a

Date: 18.04.2025

ENGLISH



Individual weight, OIML Class E2, 1 g, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-01-2020a
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Trade parts list element
EAN code (Model)	404576 100352 1
Successor Model	317-01; 317-01-2022i
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	0,03 mg
Case aluminium	317-010-600
Case plastic	317-020-400
Case wood	317-010-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

OIML class	E2
OIML conform	yes
Density (kg/m ³)	7950
Volume (cm ³)	0,125

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×6 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-331
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-010-600	Aluminium Protective Box KERN 317-010-600
317-010-100	Wooden Box for Set of Weights KERN 317-010-100
317-020-400	Plastic Weight Box KERN 317-020-400

SERVICES

Model	Description
952-331	Verification KERN 952-331
962-331R	DAkkS Recalibration KERN 962-331R
962-331	DAkkS Calibration KERN 962-331
920-33	Marking by etching KERN 920-33
920-331	Marking by etching (last digits of the serial number) KERN 920-331
921-331	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-331