

PRODUCT BROCHURE

KERN KUP-13

Date: 18.04.2025

ENGLISH



Extension box for connecting up to three interfaces in parallel

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Extension box for connecting up to three interfaces in parallel
- The KERN Universal Port (KUP) enables the exchange of data and control commands between balances and the network or PC, mobile end unit, printer, etc. using several interfaces: RS 232, USB device, Ethernet, Bluetooth and WLAN. In this way we ensure the perfect integration of your KERN balance into production or processing chains, which gives the benefit of a simple, consistent workflow
- In accordance with the „Plug and Play“ principle, the KUP for RS 232 or USB, Ethernet, Bluetooth or WLAN will only be inserted into the DB 15 socket which is integrated into the device. It is therefore possible to retrofit or replace interfaces conveniently without having to open the balance casing and no complicated installation. This option therefore offers a high level of flexibility at low additional cost

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	KUP-13
Model-Serie	KUP
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	4045761254008
Customs tariff number	90160090

Packaging component - 100 g
by weight - styrofoam

Delivery time 1 d

FUNCTIONS

Standard



Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	225×145×40 mm
Gross weight	0.50 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	100 g

ACCESSORIES

Model	Description
KUP-01	Adapter KERN KUP-01
KUP-03	Adapter KERN KUP-03
KUP-04	Adapter KERN KUP-04
KUP-05	Adapter KERN KUP-05
KUP-06	Adapter KERN KUP-06
KUP-08	Adapter KERN KUP-08