

KERN 317-04

Individual weight, OIML Class E2, 10 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	317-04
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761312524
Predecessor Model	317-04-2020a; 317-04-2022i
Successor Model	317-04-2022i
Planned start date	2022-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 g
Tolerance (OIML)	0,06 mg
Case aluminium	317-040-600
Case plastic	317-040-400
Case wood	317-040-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	10 g
OIML class	E2
OIML conform	Yes

Functions

Standard

Tolerance (OIML) 0,06 mg

Density (kg/m³) 7950

Volume (cm³) 1,25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 10×10×18 mm

Gross weight 0.05 kg

Net weight 0.05 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Technical data - Verification

Verification possible Yes

Optional Services

DAkkS Calibration 962-334

KERN 317-04

Individual weight, OIML Class E2, 10 g, Knob, stainless steel polished

