


PRODUCT BROCHURE

KERN KUM-08

Date: 18.04.2025

ENGLISH



Analogue module

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- The KERN Universal Modules (KUM) enable the exchange of data and control commands between balance and network or PC, mobile end unit, printer etc. using several interfaces: RS 232, USB device, Ethernet, Bluetooth and WLAN. In this way we ensure the perfect integration of your KERN balance into production or processing chains, which gives the benefit of a simple, consistent workflow
- The KUM interface modules for RS 232, USB, Ethernet, Bluetooth and WLAN are plugged straight onto the board to suit individual choice and this means that the IP protection of the housing remains guaranteed
- The KUM-08 analogue module converts the measured weight proportionately into an analogue signal. This signal can then be used to control an analogue system, which reacts to the weight

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	KUM-08
Model-Serie	KUM
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	404576 1476783
Customs tariff number	90160090

Packaging component - 100 g
by weight - styrofoam

Delivery time 1 d

FUNCTIONS

Standard



Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	65×60×100 mm
Gross weight	0.50 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	100 g