


PRODUCT BROCHURE

SAUTER AC 12R

Date: 20.04.2025

ENGLISH



Universal 2-jaw grip for tension and rupture tests

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Scope of delivery: 1 piece

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	AC 12R
Model-Serie	AC
Assortment	Sauter
Product Type	Force measurement
Product group	Clamps & clamping tools
Product usage type	Trading parts list element;Accessories/Services
EAN code (Model)	404576 1122727
Customs tariff number	90278990

Technical data - Primary

Max. Force	5 kN
Max. Clamp width	25 mm

Technical data - Construction

Material	steel
Dimension housing (W×D×H)	110×80×25 mm

Technical data - Measuring system

Mounting threads	External thread M10
------------------	---------------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×100×10 mm
------------------------------	---------------

Gross weight	0.50 kg
--------------	---------

Net weight	0.50 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

UN number	No dangerous goods
-----------	--------------------

Packaging component - by weight - plastics	10 g
--	------

Delivery time	1 d
---------------	-----

FUNCTIONS

Standard



ACCESSORIES

Model	Description
AFM 06	Adapter KERN AFM 06
AFM 03	Sleeve SAUTER AFM 03
AFM 16	Adapter KERN AFM 16
AFM 11	assembly and repair AFM 11