

PRODUCT BROCHURE

KERN YPS-05

Date: 18.04.2025

ENGLISH



Antivibration plate

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- KERN antivibration plates {=VARVAL("{#1}")} have been constructed to absorb vibrations and oscillations, which would otherwise distort the weighing result from highly-precision analytical and precision balances. Also ideal for vibration-free positioning of microscopes
- The antivibration plate is made from polished granite and is positioned on absorbing rubber components.
- You can work in an ergonomic way at the antivibration plate whether you are in a standing or seated position.
- Compact size, practical for small spaces
- All surfaces are hygienic and easy to clean
- Suitable for all KERN analytical balances, precision balances and microscopes with appropriate dimensions
- Overall dimensions W×D×H 565×450×60 mm
- Net weight approx. 42 kg

Ready to send in 2 days. Dispatch complete on pallet.
dimensions, gross weight, please ask for shipping costs

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YPS-05
Model-Serie	YPS
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	4045761253568
Customs tariff number	84239010

Net weight	42 kg
Shipping method	Freight
Packaging component - by weight - cardboard	4.5 kg
Packaging component - by weight - plastics	4.5 kg
Delivery time	2 d

FUNCTIONS

Standard



Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	800×600×400 mm
Gross weight	50 kg

SPARE PARTS

Model	Description
RTYPS-8001	Rubber shock absorber KERN RTYPS-8001