

PRODUCT BROCHURE

OPTICS OBB-A240 1

Date: 19.04.2025

ENGLISH



Trinocular lens barrel

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Microscope tube OBB-A2401

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A2401
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761172760
Customs tariff number	90119000

Storage Temp Max 40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×185×165 mm
Gross weight	2 kg
Net weight	1.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	600 g
Delivery time	1 d

ENGLISH

Technical data - Primary

Tube	Binocular
------	-----------

Technical data - Construction

Tube type	Siedentopf
Tube angle	30°
Eye distance (max)	75 mm
Eye distance (min)	52 mm
Adjustable dioptr	one-sided

Technical data - Environmental conditions

Storage Temp Min	-5 °C
------------------	-------