


KERN TVS 10KN100-LD

Practical set for force-displacement-measurements in laboratory and industry



Article codes, product groups & Status information

Model article-code	TVS 10KN100-LD
Model-Serie	TVS-LD
Assortment	Sauter
Product Type	Force measurement
Product group	Test stands (motorised)
Product usage type	Main article
EAN code (Model)	4045761343573
Planned start date	2023-2HJ
Customs tariff number	90249000

Technical data - Primary

Max. Force	10 kN
------------	-------

Technical data - Construction

Dimensions of housing	250×400×1550 mm
Material housing	Coated metal
Guide column(s) - distance	215 mm
Mounting threads	4x Gewindebohrung M6; Gewindebohrung M12

Technical data - Display

Display type	LED segments
--------------	--------------

Technical data - Power supply

Power supply	Power supply unit
Input voltage power supply / max. current	230 V AC, 50 Hz

Power socket for power supply unit

Cold device cable

Protection Class

Protection class I - protective conductor

Technical data - Measuring system

Traversing speed - min.	1 mm/min
Traversing speed - max.	200 mm/min
Max. Travel distance	210 mm

Technical data - Functions

Electrical limit switches	Yes
Repeat function	Yes
Speed display	Yes

Technical data - Environmental conditions

Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	520×345×1610 mm
Gross weight	105 kg
Net weight	65 kg
Shipping method	Freight
UN number	No dangerous goods
Packaging component - by weight - plastics	300 g
Packaging component - by weight - wood	14 kg

KERN TVS 10KN100-LD

Practical set for force-displacement-measurements in laboratory and industry



Packaging component - by weight - foam 1.8 kg

Packaging component - by weight - metal 2 kg

Delivery time 2 d

Functions

Standard

