

PRODUCT BROCHURE

OPTICS OBB-A1421

Date: 18.04.2025

ENGLISH



Dark field condenser 0,85-0,91

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Darkfield Condenser OBB-A1421

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1421
Model-Serie	OBB-AC
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761179301
Customs tariff number	90119000

Gross weight	0.80 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300 g
Delivery time	1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Construction

Condenser Type	Dark-field condenser
Condenser aperture	0.85-0.91
Condenser mounting	centerable

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×70 mm
------------------------------	-------------

SPARE PARTS

Model	Description
OBB-A1422	Microscope condenser KERN OBB-A1422
OBB-A1423	Darkfield condenser OBB-A1423