

OPTICS OBB-A 1579

Date: 26.02.2025

ENGLISH



Binocular lens barrel

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Microscope tube OBB-A1579

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1579
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761246478
Predecessor Model	OBB-A1548
Customs tariff number	90119000

Technical data - Primary

Tube	Trinocular
------	------------

Technical data - Construction

Tube type	Butterfly
Tube angle	30°
Tube rotation	yes
Eye distance (max)	75 mm
Eye distance (min)	50 mm
Adjustable dioptr	one-sided

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×185×165 mm
Gross weight	2 kg
Net weight	1.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	600 g
Delivery time	1 d

FUNCTIONS

Standard

