

PRODUCT BROCHURE

OPTICS OBB-A3207

Date: 18.04.2025

ENGLISH



Achromatic objective lens, 100 x / 1,25 (oil) (sprung) W.D. (0,2 mm)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Achromatic objective lens, 100 x / 1,25 (oil) (sprung) W.D. (0,2 mm) OBB-A3207

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A3207
Model-Serie	OBB-AO
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761237469
Customs tariff number	90119000

Shipping method Parcel service

Packaging component - by weight - cardboard 300 g

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Construction

Dimensions (W×D×H)	36×25×25 mm
--------------------	-------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	65×40×40 mm
Gross weight	0.60 kg
Net weight	0.30 kg