

## PRODUCT BROCHURE

---

# OPTICS OBB-A1472

Date: 19.04.2025

ENGLISH



Binocular lens barrel

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

Tube Binocular OBB-A1472

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                 |
|-----------------------|---------------------------------|
| Model article-code    | OBB-A1472                       |
| Model-Serie           | OBB-AT                          |
| Assortment            | Optics                          |
| Product Type          | Microscopes                     |
| Product group         | Transmitted light microscopes   |
| Product usage type    | Spare part;Accessories/Services |
| EAN code (Model)      | 4045761197473                   |
| Customs tariff number | 90119000                        |

### Technical data - Primary

|      |           |
|------|-----------|
| Tube | Binocular |
|------|-----------|

### Technical data - Construction

|                    |            |
|--------------------|------------|
| Tube type          | Siedentopf |
| Tube angle         | 45°        |
| Tube rotation      | yes        |
| Eye distance (max) | 75 mm      |
| Eye distance (min) | 55 mm      |
| Adjustable dioptre | both-sided |

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Storage Temp Max | 40 °C |

### Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 265×185×165 mm |
|------------------------------|----------------|

|                                               |                |
|-----------------------------------------------|----------------|
| Gross weight                                  | 1.6 kg         |
| Net weight                                    | 1.2 kg         |
| Shipping method                               | Parcel service |
| Packaging component -<br>by weight - plastics | 300 g          |
| Delivery time                                 | 1 d            |

## FUNCTIONS

Standard

