

# Software EasyTouch

## SET-51

# Audit Log

EasyTouch Audit Logs – for critical weighing processes

---

#### CONTACT



Address

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen  
Germany



Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0  
Fax : +49 - [0] 7433 - 9933 - 149



Internet

Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## DESCRIPTION

- Note: Prerequisite for this set is the basic program SET-01 Base
- The SET-51 Audit Logs function logs all actions and activities in EasyTouch to guarantee full monitoring
- Audit Trail: All actions will be stored in the system in a tamper-proof memory and are audit trail compliant, to fulfil the highest security standards. This is vital in critical sectors such as the pharmaceutical industry, the manufacture of food-related products and the automotive and aviation industries, where comprehensive documentation of weighing processes is essential
- Login and Logout: Each login and logout in EasyTouch will be logged with the user name and time stamp
- User administration: All actions relating to the creation, modification or deletion of users will be reliably logged in the Audit Log
- User roles: User role administration will be documented, including creation, modification or deletion of roles
- Master data management: Changes to master data such as creation, modification or deletion will also be recorded
- Logging of measuring data: All measurement and weighing results from all active EasyTouch functions (here from Batch and Statistics) are also recorded in audit logs
- Software settings: Modifications to the general settings, including the Admin area, will be noted in the Audit Log

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	SET-51
Model-Serie	ETS
Assortment	KERN
Product Type	Software
Product group	EasyTouch
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1486355
Customs tariff number	85234910

Delivery time 1 d

### FUNCTIONS

Standard



Option



### Technical data - Packing & Shipping