


PRODUCT BROCHURE

SAUTER AC 18

Date: 20.04.2025

ENGLISH



Screw tension clamp with 2 wide jaw grips for tension and extraction tests

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

Screw tension clamp with jaws AC 18 (2 pieces)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	AC 18
Model-Serie	AC
Assortment	Sauter
Product Type	Force measurement
Product group	Clamps & clamping tools
Product usage type	Accessories/DL
EAN code (Model)	404576 11224 13
Customs tariff number	90319000

Technical data - Primary

Max. Force	5 kN
Max. Clamp width	60 mm

Technical data - Construction

Material	steel
Dimension housing (W×D×H)	160×100×140 mm

Technical data - Measuring system

Mounting threads	External thread M10
------------------	---------------------

Technical data - Environmental conditions

Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	160×100×140 mm
Gross weight	0.80 kg
Net weight	0.75 kg

Shipping method	Parcel service
UN number	No dangerous goods
Packaging component - by weight - plastics	10 g
Delivery time	1 d

FUNCTIONS

Standard



ACCESSORIES

Model	Description
AFM 03	Sleeve SAUTER AFM 03
AFM 06	Adapter KERN AFM 06
AFM 11	assembly and repair AFM 11
AFM 16	Adapter KERN AFM 16