

PRODUCT BROCHURE

SAUTER AD 0305

Date: 20.04.2025

ENGLISH



Image not
available

Bending device only for force-measuring devices up to max. = 2,5 kN

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	AD 0305
Model-Serie	AD
Assortment	Sauter
Product Type	Force measurement
Product group	Clamps & clamping tools
Product usage type	Trading parts list element;Accessories/Se rvices
EAN code (Model)	4045761175587
Customs tariff number	90278990

UN number	No dangerous goods
-----------	--------------------

Packaging component - by weight - plastics	0 g
---	-----

Delivery time	14 d
---------------	------

ENGLISH

Technical data - Construction

Material	steel
Dimension housing (W×D×H)	200×40×30 mm

Technical data - Environmental conditions

Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	200×40×30 mm
Gross weight	1.2 kg
Net weight	1.2 kg
Shipping method	Parcel service