

PRODUCT BROCHURE

OPTICS OBB-A1341

Date: 19.04.2025

ENGLISH



Trinocular lens barrel

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Microscope tube trinocular OBB-A1341

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1341
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761168763
Customs tariff number	90119000

Technical data - Primary

Tube	Trinocular
------	------------

Technical data - Construction

Tube type	Siedentopf
Tube angle	30°
Tube rotation	yes
Eye distance (max)	75 mm
Eye distance (min)	50 mm
Adjustable dioptr	one-sided

Technical data - Environmental conditions

Storage Temp Min -5 °C

Storage Temp Max 40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 190×130×100 mm

Gross weight 1.8 kg

Net weight 1.2 kg

Shipping method Parcel service

Packaging component - by weight - plastics 600 g

Delivery time 1 d

FUNCTIONS

Standard

