

PRODUCT BROCHURE

OPTICS OBB-A1630

Date: 19.04.2025

ENGLISH



Fluorescence unit for inverted microscope

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

100W LED fluorescent unit (B / G) OBB-A1630, , , length

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A 1630
Model-Serie	OBB-AFL
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761268128
Customs tariff number	90119000

Shipping method Parcel service

Packaging component - by weight - cardboard 500 g

Packaging component - by weight - plastics 500 g

Delivery time 1 d

Technical data - Verification

WEEE yes

ENGLISH

Technical data - Power supply

Protection Class	Protection class I - protective conductor
------------------	---

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	410×390×310 mm
Gross weight	9 kg
Net weight	8 kg

SPARE PARTS

Model	Description
OBB-A1629	Central section for 4xd HBO KERN OBB-A1629