

KERN 317-01

Individual weight, OIML Class E2, 1 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	317-01
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761312432
Predecessor Model	317-01-2020a; 317-01-2022i
Planned start date	2022-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	0,03 mg
Case aluminium	317-010-600
Case plastic	317-020-400
Case wood	317-010-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	1 g
OIML class	E2
OIML conform	Yes
Tolerance (OIML)	0,03 mg

Functions

Standard

Density (kg/m³) 7950

Volume (cm³) 0,125

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×6 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible Yes

Optional Services

DAkkS Calibration 962-331

KERN 317-01

Individual weight, OIML Class E2, 1 g, Knob, stainless steel polished

