

# KERN 356-88

Individual weight, OIML Class M2, 20 kg, Block, cast iron lacquered



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 356-88                         |
| Model-Serie           | 356                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M2 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761248274                  |
| Predecessor Model     | 356-88-2020a                   |
| Planned start date    | 2020-1HJ                       |
| Customs tariff number | 84239090                       |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 20 kg       |
| Tolerance (OIML) | 3 g         |
| Case aluminium   | 346-080-600 |

## Technical data - Construction

|            |                     |
|------------|---------------------|
| Dimensions | 230×115×139 mm      |
| Material   | Cast iron lacquered |

## Technical data - Measuring system

|                              |       |
|------------------------------|-------|
| Nominal value                | 20 kg |
| OIML class                   | M2    |
| OIML conform                 | Yes   |
| Tolerance (OIML)             | 3 g   |
| Density (kg/m <sup>3</sup> ) | 7100  |

## Functions

Standard

Volume (cm<sup>3</sup>) 2816,9014

## Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 230×115×139 mm |
| Gross weight                               | 20 kg          |
| Net weight                                 | 20 kg          |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

## Technical data - Verification

|                       |     |
|-----------------------|-----|
| Verification possible | Yes |
|-----------------------|-----|

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-645 |
|-------------------|---------|

# KERN 356-88

Individual weight, OIML Class M2, 20 kg, Block, cast iron lacquered

