

PRODUCT BROCHURE

OPTICS OBB-A3205

Date: 18.04.2025

ENGLISH



Achromatic objective lens, 40 x /0,65 (sprung) W.D. (0,6 mm)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Achromatic objective lens, 40 x /0,65 (sprung) W.D. (0,6 mm) OBB-A3205

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A3205
Model-Serie	OBB-AO
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	404576 1237445
Customs tariff number	90119000

Shipping method Parcel service

Packaging component - by weight - cardboard 300 g

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Construction

Dimensions (W×D×H) 33×25×25 mm

Technical data - Environmental conditions

Storage Temp Min -5 °C

Storage Temp Max 40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 65×40×40 mm

Gross weight 0.60 kg

Net weight 0.30 kg