

## PRODUCT BROCHURE

---

# KERN SET-261

Date: 19.04.2025

ENGLISH

# Software EasyTouch SET-261 **Audio Management**

EasyTouch Audio Management - EasyTouch Audio Management – Individual voice and sound files

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- **Audio Management:** With this option, individual voice or sound files can be stored in the system for specific events. As soon as the event occurs, the system plays the individual sound file instead of the standard sound file in this case. In this way, for example, voice output can be stored for workshops for the blind, such as “too light”, “OK” or “too heavy”. Furthermore, the volume of the sound output can be adjusted individually

## SPECIFICATIONS





### Article codes, product groups & Status information

Model article-code	SET-261
Model-Serie	ETS
Assortment	KERN
Product Type	Software
Product group	EasyTouch
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761246157
Customs tariff number	85234910

### Technical data - Packing & Shipping

Delivery time	1 d
---------------	-----

## FUNCTIONS

Standard
   
Option
