

PRODUCT BROCHURE

KERN YKA-14

Date: 18.04.2025

ENGLISH



Mains adapter , standard EURO/US/AUS/UK/CH

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Power adapter YKA-14 (15 V 4 A | EURO, UK, US, CH)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKA-14
Model-Serie	YKE
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761185012
Customs tariff number	85044095

Packaging component - 23 g
by weight - plastics

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Power supply

Mains adapter/adaptor included	CH;EURO;UK;US
Input Voltage	100 V - 240 V AC, 50 / 60 Hz

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	90×70×70 mm
Gross weight	0.35 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	23 g