

PRODUCT BROCHURE

KERN YKA-11

Date: 18.04.2025

ENGLISH



Mains adapter 100 - 240 V~ 50/60 Hz 190 mA, Standard EURO/US/AUS/UK/CH

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Power adapter YKA-11 (12 V, 500 mA | 100 - 240 V~ 50/60 Hz 190 mA | EURO, AUS, UK, US, CH)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKA-11
Model-Serie	YKA
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761177291
Customs tariff number	85044095

Delivery time

1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Power supply

Mains adapter/adaptor included	AUS;CH;EURO;UK;US
Input Voltage	100-240V~ 50/60Hz 190mA

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	90×55×80 mm
Gross weight	0.20 kg
Net weight	0.15 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	36 g