

PRODUCT BROCHURE

KERN 356-89

Date: 19.04.2025

ENGLISH



Individual weight, OIML Class M2, 50 kg, Block, cast iron lacquered

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	356-89
Model-Serie	356
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761248304
Predecessor Model	356-89-2019a
Customs tariff number	84239090

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	8 g

Technical data - Construction

Dimensions (W×D×H)	310×155×192 mm
Material	cast iron lacquered

Technical data - Measuring system

Nominal value	50 kg
OIML class	M2
OIML conform	yes
Tolerance (OIML)	8 g
Density (kg/m ³)	7100
Volume (cm ³)	7042,2535

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	310×155×192 mm
Gross weight	50 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-646
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-290	Gloves KERN 317-290

SERVICES

Model	Description
962-646	DAkkS Calibration KERN 962-646
920-30	Marking with stroke number KERN 920-30
962-646R	DAkkS Recalibration KERN 962-646R
952-646	Verification KERN 952-646
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301

ENGLISH

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium
356-86	5 kg	800 mg	346-060-600
356-76	5 kg	800 mg	346-060-600
356-77	10 kg	1,6 g	346-070-600
356-87	10 kg	1,6 g	346-070-600
356-88	20 kg	3 g	346-080-600
356-78	20 kg	3 g	346-080-600
356-79	50 kg	8 g	-
356-89	50 kg	8 g	-