

PRODUCT BROCHURE

KERN YPS-03

Date: 19.04.2025

ENGLISH



Weighing table to absorb vibrations and oscillations, which would otherwise distort the weighing result

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- The KERN YPS-03 weighing table has been constructed to absorb vibrations and oscillations, which would otherwise distort the weighing result
- The YPS-03 weighing table consists of:
 - a heavy, polished granite slab (24 kg) and is positioned on absorbing rubber components and is used as a surface on which to stand the analytical balance
 - a work table, the work surface of which surrounds the slab and offers many storage options. The sturdy powder-coated steel frame has four height-adjustable feet, which can be used to accurately level the weighing table if the floor is at all uneven
- You can work in an ergonomic way at the weighing table whether you are in a standing or seated position. The inlet on the front of the table offers enough space for you to sit down and operate the scale without touching the weighing table
- All surfaces are hygienic and easy to clean
- Suitable for all KERN analytical balances with overall dimensions =W×D 400×450 mm
- Overall dimensions W×D×H 800×600×785 mm (height is adjustable)
- Dimensions of granite slab, W×D×H 410×270×70 mm
- Net weight approx. 46 kg

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YPS-03
Model-Serie	YPS
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	404576 1246898
Predecessor Model	YPS-03-2019a
Customs tariff number	84239010

Technical data - Construction

Overall dimensions mounted (W×D×H)	800×600×785 mm
------------------------------------	----------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	800×600×300 mm
Gross weight	55 kg
Net weight	46 kg
Shipping method	Freight
Packaging component - by weight - cardboard	3750 g
Packaging component - by weight - plastics	3750 g

Delivery time 2 d

FUNCTIONS

Standard



SPARE PARTS

Model	Description
RTYPS-9001	Adjustable base foot KERN RTYPS-9001
RTYPS-8002	Rubber shock absorber KERN RTYPS-8002