

PRODUCT BROCHURE

KERN MPS-A03

Date: 19.04.2025

ENGLISH



Power supply cable

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Power adapter MPS-A03 (15 V, 300 mA | 120 V | US)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	MPS-A03
Model-Serie	MPS
Assortment	KERN
Product Type	Medical products
Product group	Personal scales
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761104310
Successor Model	YKA-29-2021a
Customs tariff number	85044095

Packaging component - 8 g
by weight - plastics

Delivery time 1 d

Technical data - Verification

Verification approval yes
conforming to
2014/31/EU

ENGLISH

Technical data - Power supply

Mains adapter/adapter included	US
Input Voltage	120V

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	89×70×55 mm
Gross weight	0.30 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	8 g